

Installation of Comark USB Drivers


Please download the latest drivers from:

<http://www.comarkltd.com/software>

This procedure is valid for Windows XP and Windows 7.


Run ComarkVCPInstaller.exe from the "Comark USB Drivers-V3_00" folder, the following screen should be displayed:

NOTE: Windows 7 User Account Control will ask if the user wants to allow this program to run, select **yes**.



Click Install

Once installation is complete you will see the following screen:



Please note your computer may request a restart after installation completes.

Comark Instruments
52 Hurricane Way,
Norwich, Norfolk, NR6 6JB England
Tel: 01603 (+44 1603) 256647
Fax: 01603 (+44 1603) 256744
Email: service@comarkltd.com
Website: www.comarkltd.com

Comark Instruments
PO Box 9090, Everett,
WA 98206, USA
Tel: (503) 643 5204
Fax: (503) 644 5859
Email: sales@comarkUSA.com
Website: www.comarkUSA.com